Lab Session 12, Experiment 11: Valence, the Combining Capacity of Elements and Molecules

One molar mass of an element or compound that produces or reacts with 1.008 grams of hydrogen has a valence of I. It follows that the same quantity producing or reacting with 2.016 grams of hydrogen is valence II, and producing or reacting with 3.024 grams of hydrogen is valence III, etc.

In this experiment, the magnesium metal will be reacted with concentrated hydrochloric acid to produce hydrogen gas according to the following chemical equation:

$$Mg_{(sol)} + x HCl_{(aq)} \rightarrow MgCl_{x(aq)} + \frac{x}{2}H_{2(gas)}$$

The hydrogen gas produced will be collected, its volume measured, and the mass in grams calculated using the ideal gas law:

$$PV = nRT$$

When using this equation, the units of the gas constant (R) must be consistent with the units used for the other quantities in the equation. The Handbook of Chemistry and Physics lists the value of R as

$$R = 0.08206 L \times atm \times mol^{-1} \times K^{-1}$$


Thus, pressure (P) must be expressed in atmospheres, volume (V) in liters, and temperature (T) in Kelvin.

11A Experiment

A sample of magnesium metal with an approximate mass of 0.24 grams will be reacted with excess hydrochloric acid. The hydrogen gas produced will be collected, its volume measured, and finally its mass in grams will be calculated.

Chemicals Required:	6M hydrochloric acid Note: 12M HCl _(aq) must be diluted.	Magnesium metal, ribbon or granules
Apparatus Required:	500 mL Florence flask Glass tubing Watch glass Rubber tubing	125 mL Erlenmeyer flask Metal or plastic trough 100 mL graduated cylinder Thermometer

The apparatus should be set up as shown in Figure 12.1. Make sure the bottom of the thistle tube almost touches the bottom of the Erlenmeyer flask. You may be instructed to use the end sink rather than a metal or plastic trough. Be sure that all fittings are tight. This includes the fittings of the rubber stopper in the Erlenmeyer flask, the glass tube through the stopper, and the rubber tubing to the glass tube. (The latter may require a hose clamp or wire winding.)



- 1. Record the exact mass of the magnesium as entry (a) in the data table.
- 2. Remove the two-hole stopper from the Erlenmeyer flask and add enough water the cover the bottom of the thistle tube. If the magnesium is ribbon, form it into a ball or tight coil, and place it in the flask. Reinsert the rubber stopper tightly.
- 3. Add 10 mL of 6M hydrochloric acid solution by pouring it into the thistle tube. Put a rubber stopper in the top of the thistle tube to avoid back-splashing acid or escape of hydrogen gas.
- 4. After the reaction is complete, allow a few minutes for the flask to cool. Be sure to keep the end of the rubber tubing above the water level in the Florence flask. The gas collected (H₂ and water vapor) must be at the same pressure as the barometric pressure. To equalize pressures, remove the gas delivery tube and raise or lower the Florence flask so that the water levels in the trough (or end sink) and in the Florence flask are equal. (You may have to tilt the flask to equalize levels. Be careful that the mouth of the Florence flask stays below the water level in the trough.) Cover the Florence flask opening with a watch glass, invert the Florence flask and place it on the desktop.
- 5. Using a graduated cylinder, measure the volume of water that must be added to the Florence flask to fill it. It will be assumed that this is the volume of the gas collected. It should be recorded as item (e) in the data table.
- 6. Measure the temperature of the water in the trough with your thermometer and record it as item (b) below. It will be assumed that this is the temperature of the collected gas.
- 7. Record the barometric pressure as item (c) below. This is the total pressure of the hydrogen and water vapor.

(a) Mass of magnesium	g
(b) Temperature of water in trough	°C
(c) Barometric pressure	Torr
(d) Vapor pressure of water at the temperature recorded in (b)	Torr
(e) Volume of gas collected = Volume of H_2O added	mL

11B Exercise: Converting Data to the Correct Units

1.	Calculate the Pressure	P =atm
	Pressure [atm] = Barometric Pressure	 Vapor Pressure of H₂O
	760	Torr/atm
2.	Calculate the Volume	V =L
	Volume in L = Vol	ume in mL \div 1000
3.	Calculate the Temperature	T =K
	Absolute temperature in Kelv	$ins = {}^{\circ}C + 273.15$
4.	Calculate the mass of hydrogen	Mass =g
	Calculating the grams of hydrogen produ	ced using the ideal gas law:
	$PV = Mass of hydrogen produced \div molar mass of h$	hydrogen × 0.08206 (R) × T

Rearranged:

Mass in grams of hydrogen produced = $P \times V \times Molar Mass of H_2$ 0.08206 (R) × T

5. From the generalized balanced equation: $Mg + x HCl \rightarrow MgCl_x + \frac{1}{2} x H_2$, we can calculate how many grams of hydrogen should be produced from the number of grams of magnesium reacted, as follows:

$$g_{\text{hydrogen}} = (g_{\text{Mg}}) \left(\frac{1 \text{ mole Mg}}{24.305 \text{ g Mg}}\right) \left(\frac{(x/2) \text{ moles H}_2}{1 \text{ mole Mg}}\right) \left(\frac{2.016 \text{ g H}_2}{1 \text{ mole H}_2}\right) = \frac{x}{24.11}.$$

Rearrangement yields the following ratio.

$$\frac{g_{\text{hydrogen}}}{g_{\text{Mg}}} = \frac{x}{24.11} \text{ Thus, } x = \left(\frac{g_{\text{hydrogen}}}{g_{\text{Mg}}}\right)(24.11) \text{ } x = \left(\frac{g_{\text{H}_2}}{g_{\text{Mg}}}\right)(24.11) = \underline{\qquad},$$

$$\% \text{ error} = \frac{|x-2|}{2} \text{ } X \text{ } 100 = \underline{\qquad}\%$$

Give a balanced equation for the reaction of magnesium with hydrochloric acid:

-	g Mg	g 1 mole N 24.3 g M	0		moles HCL needed
-	10 mL	10 ⁻³ L 1 mL	6.0 moles HCl L	- =	moles HCL added

6. Verification that HCl is the excess reagent by completing the following calculations:

_____ moles HCl added – _____ moles HCl needed = _____ moles of HCl unreacted

11C Alternate Calculations

1.	Moles of H ₂ = \underline{PV}_{RT} = $\underline{()()}_{PT}$ = \underline{mol}_{RT} mol.
2.	Moles of Mg = $\underline{g_{Mg}}$ = $\underline{()}$ = \underline{mol} .
3.	Mole ratio of H ₂ to Mg = $\frac{\text{Moles of H}_2}{\text{Moles of Mg}}$ = $\frac{()}{()}$ = $\frac{()}{()}$
4.	$x_{exp} = 2 \times ratio = 2 \times (___) = ___$
5.	% error = $\frac{x_{exp} - 2.00}{2.00}$ = $(\underline{) - 2.00}$ =%.

11D Exercise: Nomenclature

Give the formula of each of the following.

calcium sulfide	Iron (II) chloride
sodium sulfate	nickel (II) nitrate
calcium phosphate	magnesium hydroxide
lithium oxide	Tin (IV) bromide

Report Form 11: The Combining Capacity of Elements and Molecules

Name Partner_____Section #_____

11A Experiment

(a) Mass of magnesium	g
(b) Temperature of water in trough	° C
(c) Barometric pressure	Torr
(d) Vapor pressure of water at the temperature recorded in (b)	Torr
(e) Volume of gas collected = Volume of H_2O added	mL

11B Exercise: Converting Data to the Correct Units

P = atm 1. Calculate the Pressure V = _____L 2. Calculate the Volume T = K 3. Calculate the Temperature 4. Calculate the Mass of hydrogen Mass = _____g

5.
$$x = \left(\frac{g_{H_2}}{g_{Mg}} \right) (24.11) = \frac{6}{3} . \% \text{ error } = \frac{6}{3} .$$

Give a balanced equation for the reaction of magnesium with hydrochloric acid:

6. Verification that HCl is the excess reagent by completing the following calculations:

moles HCl added – moles HCl needed = moles of HCl unreacted

11C Alternate Calculations

- 1. Moles of $H_2 = _____mol.$ 4. $x_{exp} = 2 \times ratio = _____.$
- 2. Moles of Mg = _____mol. 5. % error = ____%.
- 3. Mole ratio of H_2 to Mg =_____.

11D Exercise: Nomenclature

calcium sulfide	Iron (II) chloride
sodium sulfate	nickel (II) nitrate
calcium phosphate	magnesium hydroxide
lithium oxide	Tin (IV) bromide